[Manualslib.com](https://manualslib.com) - The Global Manuals Library

[Manuals](#) / [Brands](#) / [Siemens Healthineers Manuals](#) / [Medical Equipment](#) / [ACUSON Origin](#) / [User's manual supplement](#) / PDF

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ACUSON Origin ACUSON Origin ICE Diagnostic Ultrasound System Supplement to the User and Reference Manuals Errata to the User and Reference Manuals

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An updated User and Reference Manual may be available for your ultrasound system. Access the most recent manual from the Siemens Healthineers Document Library at doclib.siemens-healthineers.com

ACUSON Origin

Product Version 1.1

Software Version VB20

ACUSON Origin ICE

Product Version 1.0

Software Version VB20

©2024 Siemens Medical Solutions USA, Inc.

All Rights Reserved.

Date of first issue: 2024-05

Date of revision: 2024-10

All trademarks are property of their respective owners. The following trademarks of Siemens Healthineers AG or its affiliates are registered or pending with the U.S. Patent and Trademark Office:

ACUSON, ACUSON Origin

The user documentation includes references to third-party products for informational purposes only. These references do not represent an endorsement of the third-party products referenced in the user documentation or an assumption of responsibility for the performance of third-party products.

Siemens Medical Solutions USA, Inc. reserves the right to change its products and services at any time. In addition, this publication is subject to change without notice.

Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
22010 S.E. 51st Street
Issaquah, WA 98029
U.S.A.
Phone: +1-888-826-9702
siemens-healthineers.com

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim
Germany
Phone: +49 9191 18-0
siemens-healthineers.com

Supplement to the User and Reference Manuals (English)

The supplement provides updated information for your existing user and reference manuals.

Refer to your ultrasound system's user and reference manuals for a complete set of operating instructions.

Touch Screen Overview

(New content: Page 1-16, Instructions for Use)

- 4 Additional selections based on the workflow selection and operating mode
 - During system configuration, includes selections to change the brightness of the control panel, touch screen, or keyboard without a lock key

Special Keys on the Retractable Keyboard

(New content: Page 1-24, Instructions for Use)

Key Name	Description
 (Lock)	(Available only for keyboards with a lock key) Locks or unlocks the keys on the keyboard.
 (Backlight)	(Available only for keyboards with a lock key) Adjusts the brightness of the backlighting of the keys on the keyboard.

Locking the Retractable Keyboard

(New content: Page 4-5, Instructions for Use)

(Available only for keyboards with a lock key)

To lock or unlock the keyboard:

- Press and hold the lock key for approximately six seconds.

When the keyboard is locked, the backlight of the keyboard slowly flashes.

Adjusting the Touch Screen, Control Panel, and Retractable Keyboard Lighting

(Updated content: Page 4-17, Instructions for Use)

You can adjust the brightness of the touch screen, backlighting of the controls on the control panel, and the backlighting of the keys on the retractable keyboard.

Note: The ultrasound system manufacturer assembles multiple configurations. The procedure to adjust the brightness settings is dependent on the configuration of your ultrasound system.

To adjust the brightness settings:

1. Tap the configuration button.

2. Tap **Control Panel Settings**.

3. To adjust the brightness of the control panel, rotate the required control.

Selection	Description
Touch Screen	Adjusts the brightness of the touch screen.
Orange Light	Adjusts the brightness of the lighting of an active control.
White Light	Adjusts the brightness of the white lighting of an available control.
Keyboard	(Available only for keyboards without a lock key) Adjusts the brightness of the backlighting of the retractable keyboard.

To adjust the brightness of the retractable keyboard:

- For keyboards with a lock key, press the backlight key to cycle through the settings: bright, dim, off.

- For keyboards without a lock key:
 - To increase the brightness, press and hold **Shift** and then press the up arrow on the keyboard.
 - To decrease the brightness, press and hold **Shift** and then press the down arrow on the keyboard.

Errata to the User and Reference Manuals (English)

The errata provides corrected information for the user and reference manuals.

Refer to your ultrasound system's user and reference manuals for a complete set of operating instructions.

AI Measure Option

(Updated content: Page A-21, Instructions for Use)

- Compatible transducers: 5V1, 8V3, 10V4, 5Z1, Z6T, CW2
Supported image presets: Cardiac

Environmental Requirements

Transducers

(Updated content: Page A-47, Instructions for Use)

	During Operation	During Storage or Transportation
Temperature:	---	---
5Z1, Z6T	+10°C to +35°C	-10°C to +50°C

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